

## ASRC Liquid Nitrogen Fill station Standard Operating Procedure

### Upon Entering the Fill Station Area:

- Check O<sub>2</sub> sensors are not alarming.
  - Do not enter the room if the O<sub>2</sub> sensors are alarming
- Survey station for any damage.
  - Do not use the fill station if damage is observed.
  - Report any damage to Tom Dickson ([tdickson@gc.cuny.edu](mailto:tdickson@gc.cuny.edu))

### Required items:

- PPE:
  - Insulated cryogen safety gloves
  - Eye & Face protection (Safety goggle/glasses and Face Shield)
  - Lab Coat and/or Cryogen safety apron
- Liquid nitrogen tank
- Large adjustable crescent wrench
- Personal badger account login.
  - You must be a qualified user for the liquid nitrogen system.

### Safety Considerations


- Wear required PPE during all fills.
  - It is also recommended that you wear all PPE during transportation of tanks
- Wear lab appropriate clothing and footwear
  - Closed-toe shoes, preferably of non-absorbent fabric
  - Long pants and/or skirt
- **Do not enter the room if any of the O<sub>2</sub> sensors are alarming.**
- Tank vent valve should always be pointed away from you, including during transportation.

### IN CASE OF EMERGENCY:


- **Close nitrogen supply valve** to the station.
- Contact Tom Dickson (email: [tdickson@gc.cuny.edu](mailto:tdickson@gc.cuny.edu), cell: (917) 414-4608) immediately.

# Standard Operating Procedure

1. **Enable equipment in Badger**
  - a. Open badger on the computer desktop
  - b. Sign-in using your credentials.
  - c. Select 'N2' tab on the left hand side.
  - d. Select 'Liquid Nitrogen Fill Station' to highlight blue.
  - e. Right click 'Liquid Nitrogen Fill Station'
  - f. Select 'Enable'



Note: to obtain an account and/or access to the fill station, contact Meghan Hughes ([mhughes1@gc.cuny.edu](mailto:mhughes1@gc.cuny.edu)). You must complete trainings to become a qualified user.



2. Wait for the system to initialize.
  - a. The reading on the display panel will stabilize and show a constant reading.
3. **Check the system is set to lbs.**
  - a. Check the front of the display panel. Lbs. should be highlighted.
  - b. Verify by standing on the scale. Weight will be displayed on the display panel.
4. Remove all items/debris from the scales.
  - a. This includes any hoses which may leaning on/adding weight to the scale.
5. **Ensure the scale reads zero lbs.**
6. **Roll the tank up the ramp and onto the scale.**
  - a. The vent valve should be on our left, and the liquid valve on your right.
  - b. Note: tanks typically have labelling to indicate which valve is which. If you are unsure, please contact your supplier for assistance. (or a POC?)
7. **Connect the tank valves to the fill station lines**
  - a. Connect the fill station's vent line to the vent valve of the tank first.

- b. Connect the fill line to the liquid valve of the tank.
- c. Verify that vent lines are connected to vent valve - NOT GAS or any other valve.
- 8. **Tighten the connections between the tank and fill station.**
  - a. Tighten by hand, then by wrench
  - b. Do not over-tighten the connectors, since they will freeze up during the filling step. A ¼-turn past hand tighten is normally sufficient.
- 9. **Open up the tank vent and liquid valves fully.**
  - a. Open vent valve first, then the liquid valve
- 10. Do not step on or add another weight to the scale from this point until the fill is complete
- 11. The fill station's display will show the weight of the empty tank (tare weight) plus whatever liquid nitrogen is in the tank.
- 12. **Program the fill station**
  - a. Press the "Start" button.
  - b. Enter "2"
  - c. Enter the tare weight in lbs. on the fill station
    - i. Note: tare weight is typically labelled on the tank you are using. If listed in kgs on the tank, please calculate the weight in lbs. (1kg = 2.205 lbs).
    - ii. If you are unsure of the tare weight, contact your supplier for confirmation.
  - d. Press "Enter".
  - e. Enter the fill weight in lbs.
    - i. Note: fill weight is the weight of liquid nitrogen you would like to add to the tank.
    - ii. If you are unsure of the fill weight you require, contact your supplier for verification.
    - iii. Reference for maximum fill weights:

Tank Size (L)	Suggested fill weight (lbs.)	Max. fill weight (lbs.)
240	340	425
230	325	407

- f. Press 'Enter'
- g. You will hear a 'snap' from the system as the solenoid valve opens. This indicates the start of the nitrogen fill.
- 13. If any leak from fill line connection, tighten connector until the leak stops.
  - a. If leak does not stop, hit 'stop' on the fill station panel and contact; Tom Dickson ([tdickson@gc.cuny.edu](mailto:tdickson@gc.cuny.edu)), Brian Giebel ([bgiebel@gc.cuny.edu](mailto:bgiebel@gc.cuny.edu)), Milan Begliarbekov ([mbegliarbekov@gc.cuny.edu](mailto:mbegliarbekov@gc.cuny.edu)), and/or James Aramini ([jaramini@gc.cuny.edu](mailto:jaramini@gc.cuny.edu)).
- 14. **DO NOT LEAVE THE STATION UNATTENDED DURING FILL.**
  - a. The fill station display panel will indicate the weight of liquid nitrogen added to the tank.
- 15. The fill station will automatically stop fill when the fill weight is reached.
  - a. A 'snap' from the solenoid closing indicates the end of the fill.
  - b. The scale will now show the total weight of the tank (tare weight + liquid nitrogen weight)
- 16. Close liquid valve on the tank.
- 17. Close vent valve on the tank.
- 18. **Detach hoses from the tank**
  - a. Use wrench to remove the lines from the tank.
  - b. If frozen, use a heat gun to gradually warm the connection.
- 19. **Remove the tank from the scale**
- 20. **Disable the equipment in Badger.**
  - a. Right click 'Liquid Nitrogen Fill Station'
  - b. Select 'disable'
- 21. **Close Badger** program on the PC to log out.